


Product Datasheet

BrdU antibody (orb10204)


Description	BrdU antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Brdu
Target	BrdU(Proliferation Marker)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	0.3071 kDa
CAS Number	59-14-3
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P of mice ovarian tissue (BrdU antibody)



Immunohistochemical analysis of paraffin-embedded tissue (BrdU antibody)



Immunohistochemical analysis of formalin-fixed tissue (BrdU antibody)