

Product Datasheet

Anti-phospho-GLUR4 (Ser862)/Gold [Out of stock] (orb1020823)

Description	Anti-phospho-GLUR4 (Ser862)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated synthesised phosphopeptide derived from human GluA4 around the phosphorylation site of Ser862
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	98kDa
Uniprot ID	P19493, P48058, Q9Z2W8
Expiration Date	12 months from date of receipt.