

Product Datasheet

Anti-phospho-FANCD2(Ser1404)/Gold [Out of stock] (orb1021312)

Description	Anti-phospho-FANCD2(Ser1404)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FANCD2 around the phosphorylation site of Ser1404
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	162kDa
Uniprot ID	Q80V62, Q6IV68, Q9BXW9
Expiration Date	12 months from date of receipt.