

Product Datasheet

Anti-Complement C3d fragment/Gold [Out of stock] (orb1021339)

Description	Rabbit polyclonal antibody to Anti-Complement C3d fragment/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated synthetic peptide derived from human Complement C3d fragment
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33/109/181kDa
Uniprot ID	P01024
Expiration Date	12 months from date of receipt.