


Product Datasheet

CASPR antibody (orb1023)

Description	CASPR antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human CASPR/Neurexin4 (151-250/1384 aa)
Target	CASPR
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	154 kDa
Uniprot ID	P78357
Dilution Range	ICC=1:100-500, Flow-Cyt=3ug/Test, WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...