

Product Datasheet

Anti-phospho-SMC1 (Ser966)/Gold [Out of stock] (orb1024170)

Description	Rabbit polyclonal antibody to Anti-phospho-SMC1 (Ser966)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated synthesised phosphopeptide derived from human SMC1 around the phosphorylation site of Ser966.
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	136kDa
Uniprot ID	Q9CU62, Q8NDV3, Q14683, Q920F6
Expiration Date	12 months from date of receipt.