

Product Datasheet

Anti-phospho-p38 MAPK (Tyr323)/Gold [Out of stock] (orb1025030)

Description

Anti-phospho-p38 MAPK

Species/Host

Rabbit

Conjugation

Gold

Tested

GICA, IEM, IGS

Applications
Immunogen

KLH conjugated synthetic peptide derived from human phospho-p38 MAPK (Tyr323)

Preservatives

Constituents: 0.02M TBS, pH 8.2 with 10mg/ml rAlbumin and 0.05% NaN₃, 50% glycerol. Or Lyophilized. Buffer = 0.02M TBS, pH 8.2 with 10mg/ml rAlbumin and 0.05% NaN₃. Reconstitute with sterile distilled water.

Form/Appearance

Lyophilized or Liquid

Concentration

0.4mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

41kDa

Uniprot ID

[Q16539](#), [P70618](#), [P47811](#), [O02812](#)

Dilution Range

IEM=1:20-200 IGS=1:20-200
GICA=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.