

Product Datasheet

Anti-SCARB1/Scavenger Receptor BI/Gold [Out of stock] (orb1025178)

Description

Anti-SCARB1/Scavenger Receptor

Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated synthetic peptide derived from human SCARB1/Scavenger Receptor BI
Preservatives	Constituents: 0.02M TBS, pH 8.2 with 10mg/ml rAlbumin and 0.05% NaN ₃ , 50% glycerol. Or Lyophilized. Buffer = 0.02M TBS, pH 8.2 with 10mg/ml rAlbumin and 0.05% NaN ₃ . Reconstitute with sterile distilled water.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61kDa
Uniprot ID	P97943, Q61009, Q8WTV0
Expiration Date	12 months from date of receipt.