

Product Datasheet

Anti-phospho-C-Myc(Ser62)/Gold [Out of stock] (orb1025427)

Description	Anti-phospho-C-Myc(Ser62)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human C-Myc around the phosphorylation site of Ser62
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49kDa
Uniprot ID	P01106, P01108, P09416
Expiration Date	12 months from date of receipt.