

Product Datasheet

Anti-phospho-MLKL (Ser358)/Gold [Out of stock] (orb1025509)

Description	Mouse monoclonal antibody to Anti-phospho-MLKL (Ser358)/Gold
Species/Host	Mouse
Conjugation	Gold
Immunogen	KLH conjugated synthetic peptide derived from human phospho-MLKL (Ser358)
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	54kDa
Uniprot ID	Q9D2Y4, Q8NB16
Expiration Date	12 months from date of receipt.