

Product Datasheet

Anti-phospho-PDCD4 (Ser67)/Gold (orb1026390)

Description

Anti-phospho-PDCD4 (Ser67)/Gold

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Gold

Tested

ICC

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human PDCD4 around the phosphorylation site of Ser67

Target

PDCD4

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

0.4mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

51kDa

Uniprot ID
Q53EL6
Dilution Range

IEM=1:20-200 IGS=1:20-200 GICA=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.