

Product Datasheet

Anti-Histone H4 (Acetyl K5)/Gold (orb1026772)

Description	Anti-Histone H4 (Acetyl K5)/Gold
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Gold
Tested Applications	ICC
Immunogen	KLH conjugated synthesised acetylpeptide derived from human Histone H4 around the acetylation site of Acetyl Lys5
Target	Histone H4 (Acetyl K5)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11kDa
Uniprot ID	P62805
Expiration Date	12 months from date of receipt.