

Product Datasheet

Anti-phospho-NMDAR2B (Tyr1472)/Gold [Out of stock] (orb1026844)

Description

Anti-phospho-NMDAR2B

Species/Host	Rabbit
Conjugation	Gold
Tested Applications	GICA, IEM, IGS
Immunogen	KLH conjugated Synthesised phosphopeptide derived from Rat NMDAR2B around the phosphorylation site of Tyr1472
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	180kDa
Uniprot ID	Q00960, Q01097, Q13224
Dilution Range	IEM=1:20-200 IGS=1:20-200 GICA=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.