

Product Datasheet

Anti-Dopamine Receptor D1/Gold [Out of stock] (orb1027828)

Description	Anti-Dopamine Receptor D1/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated synthetic peptide derived from human Dopamine Receptor D1
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Uniprot ID	P18901, Q61616, P21728
Dilution Range	IEM=1:20-200 ICA=1:20-200 ChIP=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.