



Product Datasheet

CD105 antibody (orb10285)


Description	CD105 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD105. Please contact us for the exact immunogen sequence. The peptide is available as orb374706.
Target	CD105
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
MW	70 kDa
Uniprot ID	P17813, P37176, Q63961
NCBI	000118, 31, 000109, 11
Entrez	2022
Dilution Range	WB: 1:100-600, IHC-P: 1:200 (based on 0.5mg/ml)
Expiration Date	12 months from date of receipt.



IHC-P staining of human breast cancer ti...



IHC-P image of human breast cancer tissu...



Immunohistochemical staining of pig hear...