



Product Datasheet

Sialoadhesin antibody (orb10301)


Description	Sialoadhesin antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Sialoadhesin. Please contact us for the exact immunogen sequence. The peptide is available as orb374709.
Target	Sialoadhesin
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	183 kDa
Uniprot ID	Q62230
Entrez	6614
Dilution Range	IHC-P: 1:50-300, WB: 1:100-2000
Expiration Date	12 months from date of receipt.



IHC-P staining of rat skin tissue using ...



Immunohistochemical staining of paraffin...



Western blot analysis of rat testis (lan...