



Product Datasheet

CD20 antibody (orb10307)


Description	CD20 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD20. Please contact us for the exact immunogen sequence. The peptide is available as orb374711.
Target	CD20
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	33 kDa
Uniprot ID	P19437, P11836
NCBI	21, 31, 021950, 068769
Entrez	931
Dilution Range	WB: 1:100-500, IHC-P: 1:100-800 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Western blot analysis of human thyroid t...



IHC-P image of human endometrial cancer ...



WB analysis of mouse lymph node (lane 1)...