



Product Datasheet

CD24 antibody (orb10309)


Description	CD24 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD24 Propeptide (60-80/80 aa)
Target	CD24
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8 kDa
Uniprot ID	P24807
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



Western blot analysis of human colon car...



rat heart lysates(Lane2) using CD24 anti...



Immunohistochemical analysis of formalin...