



Product Datasheet

CD31 antibody (orb10314)


Description	CD31 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD31. Please contact us for the exact immunogen sequence. The peptide is available as orb374712.
Target	CD31
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	82 kDa
Uniprot ID	Q3SWT0 , Q08481 , P16284
NCBI	000433 , 41 , 000442
Entrez	5175
Dilution Range	WB:2 µg/ml, IHC-P:1:500, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



IHC-P image of mouse kidney tissue using...



Immunohistochemical staining of paraffin...



IHC-P staining of mouse lung tissue using...