

Product Datasheet

CD38 antibody (orb10319)

Description

CD38 antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human CD38. Please contact us for the exact immunogen sequence. The peptide is available as orb12885.

Target

CD38

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

34 kDa

Uniprot ID

[P56528](#), [Q64244](#)

NCBI

[21](#), [031672](#), [007646](#), [41](#)

Entrez



[12494](#)

Dilution Range


IHC-P: 1:100-600, IF/ICC: 1:100-600, WB: 1:200-1500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunofluorescence analysis of guinea pi...



IHC-P staining of guinea pig pancreas ti...