

Product Datasheet

**Anti-phospho-P70 S6 Kinase beta 2 (Ser370)/Gold
[Out of stock] (orb1031924)**

Description

Anti-phospho-P70 S6 Kinase beta 2

Species/Host

Rabbit

Conjugation

Gold

Tested

GICA, IEM, IGS

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human P70 S6 Kinase beta 2 around the phosphorylation site of Ser370

Preservatives

0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

0.4mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

70kDa

Uniprot ID

Q9Z1M4, Q9UBS0

Dilution Range

IEM=1:20-200 IGS=1:20-200
GICA=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.