



Product Datasheet

PCSK9 antibody (FITC) (orb103203)

Description	PCSK9 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human PCSK9/NARC1
Target	NARC1/PCSK9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NBP7
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat Brain



Immunohistochemical staining of Mouse Brain