

## Product Datasheet

**Anti-phospho-PKC delta (Ser645)/Gold [Out of stock] (orb1032477)**

**Description**

Anti-phospho-PKC delta (Ser645)/Gold

**Species/Host**

Rabbit

**Conjugation**

Gold

**Tested**

GICA, IEM, IGS

**Applications**
**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human PKC delta around the phosphorylation site of Ser645

**Preservatives**

0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

0.4mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

77kDa

**Uniprot ID**
**P28867, Q05655, P09215**
**Dilution Range**

IEM=1:20-200 IGS=1:20-200  
GICA=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

**Expiration Date**

12 months from date of receipt.