

Product Datasheet

CD8 antibody (orb10325)

Description

CD8 antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human CD8. Please contact us for the exact immunogen sequence. The peptide is available as orb374713.

Immunohistochemical staining of paraffin...

Target

CD8

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Immunohistochemical staining of human en...

Note

For research use only

Application notes

This antibody is available without BSA and NaN₃, please contact us for details

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

IHC-P image of rat spleen tissue using a...

MW

26 kDa

Uniprot ID
P10300, P01732, P05541
NCBI
21, 001081110, 001074579, 11
Entrez
12525
Dilution Range

WB: 1:200-800, IHC-P: 1:100-300 (based on 0.5 mg/ml)

Expiration Date

12 months from date of receipt.