

Product Datasheet

CD45 antibody (orb10328)

Description

CD45 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human CD45. Please contact us for the exact immunogen sequence. The peptide is available as orb374714.

Target

CD45

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

143 kDa

Uniprot ID

P08575, P06800, P04157

Entrez

5788

Dilution Range

IF/ICC: 1:100-500, WB: 1:200-1000, IHC-P: 1:100-500 (based on 0.5 mg/ml)

Expiration Date

12 months from date of receipt.