



Product Datasheet

E-selectin antibody (orb10342)

Description	E-selectin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse E-Selectin (501-612/612 aa)
Target	E-Selectin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65 kDa
Uniprot ID	P16581
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P of human gastric carcinoma tissue ...



Immunohistochemical analysis of formalin...