


## Product Datasheet

### CDKN1C antibody (orb10364)

|                            |                                                                                                    |
|----------------------------|----------------------------------------------------------------------------------------------------|
| <b>Description</b>         | CDKN1C antibody                                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                       |
| <b>Tested Applications</b> | FC, IHC-P, WB                                                                                      |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from rat Cdkn1c (291-343/343 aa)                          |
| <b>Target</b>              | Cdkn1c                                                                                             |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                              |
| <b>Form/Appearance</b>     | Liquid                                                                                             |
| <b>Concentration</b>       | 1mg/ml                                                                                             |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                    |
| <b>Note</b>                | For research use only                                                                              |
| <b>Isotype</b>             | IgG                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                         |
| <b>MW</b>                  | 35 kDa                                                                                             |
| <b>Uniprot ID</b>          | <b>P49918</b>                                                                                      |
| <b>Dilution Range</b>      | IHC-P=1:100-500, Flow-Cyt=0.2ug/test, IHC-F=1:100-500, ICC=1:100, IF=1:100-500, ELISA=1:5000-10000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                    |



Flow  
cytometry of  
SH-  
SY5Y.Primary  
Antibo...