


Product Datasheet

CEA antibody (orb10368)

Description	CEA antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CEA/CD66e/CEACAM5 (full length)
Target	CEA(B5)
Preservatives	0.01M PBS(pH7.4) with 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	150-200 kDa
Uniprot ID	P06731
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...