

Product Datasheet

Anti-Histone H4 (Acetyl K16)/Cy5.5 (orb1036884)

Description

Anti-Histone H4 (Acetyl K16)/Cy5.5

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Cy5.5

Tested

FC, IF

Applications
Immunogen

KLKLH conjugated Synthesised
acetylpeptide derived from human
Histone H4 around the acetylation site of
K16

Target

H4C1

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin,
0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid
repeated freeze/thaw cycles. The
lyophilized antibody is stable at room
temperature for at least one month and
for greater than a year when kept at -
20°C. When reconstituted in sterile pH
7.4 0.01M PBS or diluent of antibody the
antibody is stable for at least two weeks
at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

11kDa

Uniprot ID
P62805
Dilution Range

Flow-Cyt=1ug/TEST,IF=1:50-200

Expiration Date

12 months from date of receipt.