



Product Datasheet

HER2 antibody (orb10372)


Description	HER2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human HER2 (1201-1255/1255 aa)
Target	ERBB2/HER2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	185 kDa
Uniprot ID	P04626
Dilution Range	IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human ovarian tissue (ErbB2) an...



Paraffin-embedded human cervical carcinoma...