



Product Datasheet

Chromogranin A antibody (orb10377)

Description	Chromogranin A antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CHGA (421-470/457 aa)
Target	CHGA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	P26339
Dilution Range	IHC-P=1:100-500, ICC=1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of mice liver tiss...



Immunofluorescence analysis of mice panc...