


Product Datasheet

pan Cytokeratin antibody (orb10399)

Description	pan Cytokeratin antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human cytokeratins
Target	Pan Cytokeratin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42-64 kDa
Uniprot ID	P13646
Dilution Range	WB=1:500-2000, IHC-P=1:100-2000, ICC=1:100, IF=1:100-500, Flow-Cyt=1µg /test, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...