


Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody (orb1040884)

Description	SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	> 1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	46 kDa
Uniprot ID	P0DTC2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis using Rabbit Anti-...