



Product Datasheet

APG4B(Acetyl K39) antibody (orb1040885)


Description	APG4B(Acetyl K39) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human APG4B(Acetyl K39)
Target	APG4B(Acetyl K39)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	Q9Y4P1
Dilution Range	IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat brain
using
APG4B ...



IHC-P
analysis of
mouse brain
using
APG4...



IHC-P
analysis of
human brain
using
APG4...