



Product Datasheet

CLCA4 antibody (orb10419)


Description	CLCA4 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse chloride channel calcium activated 3B precursor (651-750/1044 aa)
Target	calcium-activated chloride channel CLCA4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110 kDa
Uniprot ID	Q14CN2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of rat intestinal mucosa using CLC...



Western blot analysis of mouse liver tis...