



Product Datasheet

Collagen III antibody (orb10437)


Description	Collagen III antibody
Species/Host	Rabbit
Reactivity	Canine, Human, Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Collagen alpha 1(III) chain (1301-1400/1466 aa)
Target	Collagen III
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	117 kDa
Uniprot ID	P02461
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of rat liver tissue (Collagen III ...



Immunohistochemical analysis of paraffin...