

Product Datasheet

Anti-phospho-PIK3C3(Ser164)/Cy5.5 (orb1044107)

Description

Anti-phospho-PIK3C3(Ser164)/Cy5.5

Species/Host

Rabbit

Reactivity

Human

Conjugation

Cy5.5

Tested

FC, IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PIK3C3 around the phosphorylation site of Ser164

Target

PIK3C3

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

98kDa

Uniprot ID
Q8NEB9
Dilution Range

Flow-Cyt=1ug/Test,IF=1:100-500

Expiration Date

12 months from date of receipt.