

Product Datasheet

Anti-phospho-Alpha synuclein (Tyr125)/RBITC (orb1044155)

Description	Anti-phospho-Alpha synuclein (Tyr125)/RBITC
Species/Host	Rabbit
Reactivity	Human
Conjugation	RBITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SNCA around the phosphorylation site of Tyr125
Target	phospho-Alpha synuclein (Tyr125)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16143
Dilution Range	Flow-Cyt=3ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.