

## Product Datasheet

### Collagen X antibody (orb10444)

## Description

Collagen X antibody

## Species/Host

Rabbit

## Reactivity

Bovine, Canine, Gallus, Human, Mouse, Rat

## Conjugation

Unconjugated

## Tested

ELISA, IHC-P, WB

## Applications

## Immunogen

KLH conjugated synthetic peptide derived from human Collagen X. Please contact us for the exact immunogen sequence. The peptide is available as orb374717.

## Target

Collagen X

## Form/Appearance

10 mM PBS, 0.02% sodium azide

## Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

## Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

## Note

For research use only

## Isotype

IgG

## Clonality

Polyclonal

## Purity

Polyclonal antibodies are purified by peptide affinity chromatography

## MW

66 kDa

## Uniprot ID

[Q05306](#), [Q03692](#)

## NCBI

[000493](#), [000484](#)

## Entrez


[1300](#)

## Dilution Range


IHC-P: 1:300-500, WB: 1:200-2000

## Expiration Date


12 months from date of receipt.



IHC-P staining of rat liver tissue using...



Immunohistochemical staining of rat live...



IHC-P staining of rat liver tissue using...