

Product Datasheet

Anti-Phospho-NMDAR2B (Tyr1336)/Cy7 [Out of stock] (orb1044800)

Description	Anti-Phospho-NMDAR2B (Tyr1336)/Cy7
Species/Host	Rabbit
Conjugation	Cy7
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse NMDAR2B around the phosphorylation site of Tyr1336
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	164kDa
Uniprot ID	Q00960 , Q01097 , Q13224
Expiration Date	12 months from date of receipt.