

Product Datasheet

phospho-MEK1 (Ser385)/RBITC antibody (orb1044804)

Description	phospho-MEK1 (Ser385)/RBITC antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	RBITC
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ERK1 around the phosphorylation site of Ser385
Target	phospho-MEK1 (Ser385)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Uniprot ID	P36507
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.