

Product Datasheet

Anti-phospho-EIF4G1 (Ser1238)/Cy5.5 [Out of stock] (orb1044917)

Description

Anti-phospho-EIF4G1 (Ser1238)/Cy5.5

Species/Host

Rabbit

Conjugation

Cy5.5

Tested Applications

IF

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human EIF4G1 around the phosphorylation site of Ser1238

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

175kDa

Uniprot ID

Q04637, **Q4G001**, **Q6NZJ6**

Dilution Range

IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.