

## Product Datasheet

# Anti-Phospho-Histone H2A.X (Tyr143)/RBITC (orb1044983)

<b>Description</b>	Anti-Phospho-Histone H2A.X (Tyr143)/RBITC
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	RBITC
<b>Tested Applications</b>	FC, ICC, IF
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human Histone H2AX around the phosphorylation site of Tyr143
<b>Target</b>	Phospho-H2AX (Tyr143)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P16104</b>
<b>Dilution Range</b>	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
<b>Expiration Date</b>	12 months from date of receipt.