

Product Datasheet

Phycocyanin antibody (orb10452)

Description	Phycocyanin antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	c-phycocyanin purified
Target	C-Phycocyanin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	11016-15-2
Uniprot ID	P00308
Expiration Date	12 months from date of receipt.