



Product Datasheet

C-reactive protein antibody (orb10457)


Description	C-reactive protein antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombined CRP protein
Target	CRP(C11F2)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	25 kDa
Uniprot ID	P02741
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of Human kidney using CRP...



Western blot: Lane 1: U2Os (Human) Cell; La...



Western blot: C Reactive Protein(CRP); Pri...