



Product Datasheet

CRF antibody (orb10465)


Description	CRF antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CRF (145-185/196 aa)
Target	CRF
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	5/22 kDa
Uniprot ID	P06850
Dilution Range	IHC-P=1:100-500, IF=1:200-800, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of Mouse brain tissue.



Immunohistochemical analysis of Rat brain tissue.



Immunohistochemical analysis of Rat adrenal gland tissue.