

Product Datasheet

Anti-phospho-AMPK alpha-2 (Thr172)/Cy5.5 (orb1046565)

Description	Anti-phospho-AMPK alpha-2 (Thr172)/Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Sheep
Conjugation	Cy5.5
Tested Applications	IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human AMPK alpha 2 around the phosphorylation site of Thr172
Target	phospho-AMPK alpha-2 (Thr172)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54646
Expiration Date	12 months from date of receipt.