

Product Datasheet

Anti-Phospho-p70 S6 Kinase Beta (Thr228)/Cy5.5 (orb1047075)

Description

Anti-Phospho-p70 S6 Kinase Beta

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human p70 S6 Kinase Beta around the phosphorylation site of Thr228
Target	RPS6KB2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59kDa
Uniprot ID	Q9UBS0
Expiration Date	12 months from date of receipt.