

Product Datasheet

Anti-Phospho-Insulin Receptor Beta (Tyr1355)/Cy5.5 (orb1047104)

Description

Anti-Phospho-Insulin Receptor Beta

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy5.5
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor alpha around the phosphorylation site of Tyr1355
Target	Phospho-Insulin Receptor Beta (Tyr1355)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	3
Uniprot ID	P06213
Dilution Range	Flow-Cyt=3ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.