



Product Datasheet

CTBP1 antibody (orb10472)


Description	CTBP1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CTBP1 (351-440/440 aa)
Target	CTBP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	Q13363
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot of mouse liver tissue lysat...



IHC-P of rat ovary tissue (orb10472 at 1...